

INDEX

GAS ANALYZERS

AL 2021 H₂O₂

AL 4021 HCHO

AL 5001 CO FAST

AL 5002 CO ULTRA

RESEARCH

OUR PARTNERS

CONTACT




Further information

Principle of operation


The detection of peroxides is based on the liquid phase reaction of peroxides with p-hydroxyphenylacetic acid catalyzed by peroxidase. This reaction produces a fluorescent dimer, that can be excited at 326 nm (Cd-lamp) and detected between 400 and 420 nm. The technique is sensitive to all peroxides in the solution. To distinguish between H₂O₂ and organic peroxides two parallel channels are used.

In one channel (channel B) of the instrument, H₂O₂ is destroyed selectively by catalase prior to the fluorescent detection in the instrument. The amount of H₂O₂ is then given by the difference between the signals from the two channels - signal for total peroxide (channel A) minus the signal for total per-oxide without H₂O₂ (channel B) - corrected for the destruction efficiency of the catalase solution.



The above reactions are carried out in aqueous solution of peroxides and other reagents. Therefore, for the measurement of gaseous peroxides, these have to be trapped in aqueous solution first. This is achieved in a stripping coil by pumping air and a stripping solution (pH-Buffered water free of H₂O₂) continuously at known flow rates. The air and liquid streams are afterwards separated in a glass separator and the solution is then analyzed for peroxides.


In the instrument **AL2021**, capable for measuring H₂O₂ in air and water samples, the H₂O₂ mixing ratio in air is then calculated from the concentration in solution and the ratio of air and stripping solution flow rates. The coil size and the flow rates of air and stripping solution are optimized for quantitative stripping of H₂O₂. Due to lower solubilities however, the stripping efficiency for other peroxides is lower than that for H₂O₂ and probably varies between the 60% found for methylhydro peroxide and 100% for H₂O₂. As the composition of organic peroxides in air is unknown, the signal from channel A (after destroying H₂O₂ by catalase) gives only an approximate estimate of the concentration of organic peroxides. Therefore, the instrument **can not** determine exactly the amount of organic peroxides in air but a relative measure of the concentration is indicated.



Simplified block diagram of the H₂O₂ analyzer **AL2021**

Peroxides in aqueous solutions can be measured directly by each model of the instruments. In this mode, stripping is not necessary and, consequently, either zero air has to be applied to the samples inlet or the internal zero trap has to be switched on. Model **AL2021W** has been designed for measuring H₂O₂ only in water samples. Therefore no stripping coil is required. In case of using a permeation source for automatic calibration, a gas stripping system is required for transferring gaseous H₂O₂ into solution. Automatic routine measurements require regular calibrations of the instrument. As liquid phase standards of H₂O₂ are unstable and have to be prepared freshly directly before use, the instruments can be equipped with a gas phase calibration unit and a stripping coil without stabilized gas flow. All the molecules from a permeation gas H₂O₂ source are trapped in a known amount of water thus giving known liquid phase calibration standards. The technique is described in details by Lazrus et al.:


Lazrus et al., *Automated fluorometric method for hydrogen peroxide in atmospheric precipitation*, Analytical Chemistry **57** (1985) 917 and
Lazrus et al., *Automated fluorometric method for hydrogen peroxide in air*, Analytical Chemistry **58** (1986) 594, see Appendix.



Flow diagram for the determination of the H₂O₂ concentration using two separated channels.

The technique is linear up to peroxide concentrations of appr. 5×10^{-5} mol/l (eq. to appr. 2000 µg/l). For the measurement of higher concentrations the sample is prediluted inside the instrument. In this case the sensitivity range of the instrument can be expanded even to concentrations above 3000 µmol/l (eq. to appr. 100 mg/l). The dilution system is applied also if seawater has to be analyzed. Due to chemical reactions between seawater components and the chemistry used for the detection of peroxides otherwise a strong precipitation occurs inside the instrument. Dilution of the seawater helps to minimize this effect. Nevertheless it proved to be necessary to clean the instrument by flushing with 1n HCl solution on a regular time basis. This HCl cleaning step is included automatically in the instrument.

The instrument is designed to be installed in a 19" rack. In order to ensure a steady flow of the reagents, the **AL2021** should be installed above the cooling box containing the reagents.



Scheme of tubing between the **AL 2021** and the reagents

Features of the instrument

- ▶ Provides absolute concentrations for H₂O₂ and relative values for organic peroxides.
- ▶ Operates with gas phase and liquid phase samples.
- ▶ Fast and continuous monitoring of H₂O₂ concentration with time constant of 100 sec (10-90%).
- ▶ Fully automated operation using micro controller.
- ▶ Fully automated calibration by internal H₂O₂ permeation source or liquid H₂O₂ standards.
- ▶ Rugged and simple to use. RS232 interface and analog output.

The instrument **AL2021** operates with individual measuring ranges. The liquid calibration for **AL2021** can be defined by the user within the ranges given as examples below:

- ▶ Ranges for liquid phase measurement:

0 to 30 µg/l
(free programmable) 0 to 300 µg/l
0 to 3 mg/l
- ▶ Ranges for gas phase measurement:

0 to 20 ppbV
(free programmable) 0 to 200 ppbV
0 to 2 ppmV
- ▶ H₂O₂ detection limit (liquid phase): <70 ng/l
- ▶ H₂O₂ detection limit (gas phase): <50 ppt
- ▶ Equivalent to: 2×10⁻⁹ molar
- ▶ Zeroing: Internal zero trap
- ▶ Calibration: Liquid H₂O₂ standard or internal H₂O₂ permeation source (optional)
- ▶ Signal output (analogue): 0-5 V FS. Calibration automatically sets output to 4.0 V for maximum value of chosen range.
- ▶ Signal output (digital): RS232 interface
- ▶ Operation temperature: -10°C to 40 °C
- ▶ Dimensions: Width 19": 50cm×50cm×21cm (W×L×H)
- ▶ Power requirement: 110/220 VAC or 24 VDC, ~100 W

Interferences from other substances

O ₃	30 ppt H ₂ O ₂ /100 ppb O ₃ (1:3500)
NO	12 ppt H ₂ O ₂ /100 ppb NO (1:8000)

No detectable interferences from:

SO ₂	PAN	NO ₂	Glyoxal
Isobutane	1-Butane	Formaldehyde	Benzene
Toluene	Methanol	Acetone	Methylamine
Dimethylamine	n-Butane	cis-2-Butene	trans-2-Butene
Ions, such as:	J ⁻ Cl ⁻ Br ⁻	NO ₃ ⁻ PO ₄ ⁻	Benzoate

